

A1c Go

Boronate Affinity Chromatography

POCT HbA1c Analyzer




Compact



Fully Automatic



Cartridge based



Boronate Affinity



Features

- Advanced technique: Boronate Affinity Chromatography
- Fully automatic one-step operation
- Fast and accurately answer of results in minutes
- Least interference from HbA1c Variants
- No liquid Path: no carryover & minimized maintenance
- Point-of-Care analyzer is suitable for deliver timely primary care at the physician's office

Smart Look



224mm(L)*160mm(W)*236mm(H)

Simplified Operation



Click to START



Collect 5 μ L sample, insert the pipette into R1 position



Insert the cartridge into the machine, followed by automatically testing process
Result will be displayed in 5 min

Specifications

Screen	Color touch screen	Methodology	Boronate Affinity Chromatography
Assay	HbA1c, Glycated Hemoglobin	Pack size	25 kits/box
Sample volume	5 μ L	Linear range	3.5-15.0% (NGSP)
Interval	\leq 5 min	Calibration	RFID card swipe
CV	\leq 3%	Dimensions	224mm(L)*160mm(W)*236mm(H) \leq 5 Kg
Communications	LAN port, USB, Serial Port, WiFi		
consumable	Cartridge, RFID card, disposable pipette, HbA1c Control(optional)		

CONTACT

Goldsite Diagnostics Inc.
 No. 1001, Building No.1,
 GOLDSITE Mansion,
 No. 9 Rongtian Road, Kengzi,
 Pingshan District, Shenzhen,
 China, 518122

☎ 0086 755 2689 0807
 ✉ export@goldsite.com.cn
 🌐 en.goldsite.com.cn



@ Goldsite Diagnostics